

IRON CHECK
FERROUS/FERRIC
RANGE 0.0 TO 5.0 mg/L TEST STRIPS

480025
BOTTLE OF 25 TEST STRIPS

480125
30 TEST STRIPS


This patented (U.S. Patent No. 5620658) test is the sensitive and safe alternative to wet chemical testing

Inexpensive
 Reliable
 Accurate results in under three minutes.
 Pocket sized
 On site testing
 No material safety data sheets required

No mechanical reader is required.
 Reproducible results
 Portable
 Ideal for food industry
 Safe and easy to use

Available as bottle of 25 strips (480025) or 30 individual strips (480125).

IRON RANGE	480025	480125
STEPS	11	11
PPM/mg/L	0.0 TO 5.0	0.0 TO 5.0
INCREMENTS	0.0;0.02;0.05;0.1;0.2;0.3;0.5;0.75;1.0;2.0;5.0	0.0;0.02;0.05;0.1;0.2;0.3;0.5;0.75;1.0;2.0;5.0
TEST TIME	3 minutes	3 minutes



Thermo Electron Corporation Station Road, Auchtermuchty, Fife, UK
 KY14 7DP.

Tel: +44 (0)1337 828871

Fax: +44 (0)1337 828972

E-mail: sales@thermorussell.com

<http://www.thermorussell.com>