


CW711
COMBINATION pH GLASS BODIED ELECTRODE

Applications: Best for general laboratory pH measurements.

Features: Refillable reference chamber.
 Ceramic junction allowing fast response time.

pH Range	0-14	Temp' Range/°C	0-100
Slope	95 TO 102%	Reproducibility	<0.03 pH
Pressure Range	Up to 1 Bar	Electrode Resistance	100 – 500 MΩ
Long term stability (drift) 20°C in 4.00 Buffer	<0.03pH/Hr	Eo	-30 TO +30 mV
Response Time	<30 sec	Temperature Compensation	None

(max. time for a selectable pH step between 2-12 (95% of step))



Waterproof	NO	Sensor	pH bulb
Dimensions below cap	120 x 12mm	Refillable	Yes
Cable length	1M	Cable termination	BNC
Storage solution	S00V/	Body material	Glass
Reference type	Ag/AgCl	Shelf life	12 months

The CW711 is a general glass bodied pH electrode. It can be used in most everyday laboratory samples and has a refillable reference chamber. Thermo Russell have engineered the electrode to perform consistently and reliably in general laboratory applications.

Thermo Electron Station Road, Auchtermuchty, Fife, UK KY14 7DP
 Tel: +44 (0)1337 828871 Fax: +44 (0)1337 828972
 E-mail: sales.water.uk@thermo.com Website: <http://www.thermo.com>